

# CAPP

D E N M A R K

 **ecopipette**

SINGLE CHANNEL PIPETTES

## ecopipette Innovation

- THIS PIPETTE IS BUILT TO LAST. MAKING NO COMPROMISE ON QUALITY IN PRODUCT DEVELOPMENT, RAW MATERIALS USED AND MANUFACTURING PROCESSES WE DEVELOPED THE GREENEST PIPETTE AVAILABLE ON THE MARKET. THE **ECOPIPETTE** IS NOT ONLY PACKAGED USING RECYCLABLE, BIODEGRADABLE MATERIALS, BUT CONSTRUCTED FROM SOME OF THE MOST RENEWABLE DURABLE RESOURCES AVAILABLE TO ENSURE AN ADVANCED LIFE SPAN AND GREENEST AFTER SERVICE RECYCLING.
- CAPP ECOTRADE™ PROGRAM: RETURN ANY BRAND OF OLD PIPETTES IN EXCHANGE FOR DISCOUNTED NEW ECOPIPETTES™. AT CAPP APS, WE WILL ENSURE AN ENVIRONMENTALLY SAFE WASTE DISCHARGE OF ALL OLD PIPETTES. ASK YOUR FAVOURITE AUTHORIZED CAPP DISTRIBUTOR FOR MORE INFORMATION ABOUT OUR ECOTRADE™ PROGRAMME.

CAPP **ECOPIPETTE**™ UNIQUE ADVANTAGES:

- LOCKABLE VOLUME CONTROLLER REQUIRING 50% FEWER TURNS
- DISPLAY IS FULLY VISIBLE DURING HAND OPERATIONS
- REDUCED FORCE REQUIREMENT FOR ASPIRATION AND BLOW-OUT
- THE ONLY PIPETTE WITH FIXED AND VARIABLE EXCHANGEABLE COUNTER MECHANISM
- LARGE PAD FOR EASE OF TIP EJECTION
- FULLY AUTOCLAVABLE PIPETTE BODY (NO DISASSEMBLY REQUIRED)
- DOUBLE O-RING FOR OPTIMAL TIP SEAL
- COMPATIBLE WITH ALL MAJOR PIPETTE TIP BRANDS
- ROBUST DESIGN INCLUDES METAL TIP BASE
- LOWEST INCIDENCE OF BREAKAGE RESULTING FROM ACCIDENTAL DROPPING OR TIPS MOUNTING.





## Ordering Info

### ecopipette variable volume

CAT. NO.	DESCRIPTION	INACCURACY %	IMPRECISION %	COLOR
C02-1S	CAPP PIPETTE, VARIABLE VOL. 0.2-2 UL, ECOPIPETTE	12.00 / 2.00	6.00 / 1.00	
C10-1	CAPP PIPETTE, VARIABLE VOL. 0.5-10 UL, ECOPIPETTE	2.50 / 1.00	1.80 / 0.50	
C20-1	CAPP PIPETTE, VARIABLE VOL. 2-20 UL, ECOPIPETTE	5.00 / 1.00	1.50 / 0.30	
C50-1	CAPP PIPETTE, VARIABLE VOL. 5-50 UL, ECOPIPETTE	2.00 / 1.00	1.20 / 0.30	
C100-1	CAPP PIPETTE, VARIABLE VOL. 10-100 UL, ECOPIPETTE	2.50 / 0.70	1.00 / 0.20	
C200-1	CAPP PIPETTE, VARIABLE VOL. 20-200 UL, ECOPIPETTE	2.50 / 0.60	0.70 / 0.20	
C1000-1	CAPP PIPETTE, VARIABLE VOL. 100-1000 UL, ECOPIPETTE	2.00 / 0.80	0.60 / 0.20	
C5000-1	CAPP PIPETTE, VARIABLE VOL. 1-5 ML, ECOPIPETTE	1.20 / 0.60	0.30 / 0.15	

### ecopipette fixed volume (customized volumes available on request)

CAT. NO.	DESCRIPTION	INACCURACY %	IMPRECISION %
C02-1S-001	CAPP PIPETTE, FIXED VOLUME 1 UL , ECOPIPETTE	2.50	1.80
C10-1-005	CAPP PIPETTE, FIXED VOLUME 5 UL, ECOPIPETTE	1.20	0.80
C10-1-010	CAPP PIPETTE, FIXED VOLUME 10 UL, ECOPIPETTE	1.00	0.50
C50-1-020	CAPP PIPETTE, FIXED VOLUME 20 UL, ECOPIPETTE	1.00	0.30
C50-1-025	CAPP PIPETTE, FIXED VOLUME 25 UL, ECOPIPETTE	0.80	0.35
C100-1-030	CAPP PIPETTE, FIXED VOLUME 30 UL, ECOPIPETTE	1.16	0.40
C100-1-050	CAPP PIPETTE, FIXED VOLUME 50 UL, ECOPIPETTE	1.00	0.30
C200-1-100	CAPP PIPETTE, FIXED VOLUME 100 UL , ECOPIPETTE	0.70	0.20
C200-1-200	CAPP PIPETTE, FIXED VOLUME 200 UL, ECOPIPETTE	0.60	0.20
C1000-1-500	CAPP PIPETTE, FIXED VOLUME 500 UL, ECOPIPETTE	0.80	0.15
C1000-1-1000	CAPP PIPETTE, FIXED VOLUME 1000 UL, ECOPIPETTE	0.80	0.20
C5000-1A-5000	CAPP PIPETTE, FIXED VOLUME 5000 UL, COMFORT	0.60	0.15

CXXX-1-YYY "OTHER FIXED VOL. PIPETTES(NOTE: XXX=MAX PIPETTE VOLUME; YYY=DESIRED VOLUME)"

### ecopipette Starter Packs

CAT. No.	DESCRIPTION
SKE-03-XX	THREE PIPETTES FREE OF CHOICE, STAND, RACKS OF TIPS AND TIMER