

## Disposable PCR Plate 96 well plate – 150 µL



### Product description

Riplate® PCR, 96-well, 150 µl, skirted, low profile, clear

<b>Part Number</b>	43001-0305
<b>Well volume</b>	150 µl
<b>Well bottom</b>	V-Bottom
<b>Color well</b>	clear
<b>Color frame</b>	clear
<b>Sterile</b>	no
<b>Re-usable</b>	no
<b>Patient contact</b>	no

### Classification

**Markets** EU

#### Purity




Produced under cleanroom conditions  
(ISO 14644-1 class 8)

Free from DNA, DNase, RNase, and PCR inhibitors.

### Materials

<b>Well</b>	PP – Polypropylene
<b>Frame</b>	PC – Polycarbonate

Dimensions	
Plate Length	127.76 mm
Plate Height	15.9 mm
Plate Width	85.47 mm



Packaging	
Box	50 plates
Pallet	84 boxes (4,200 plates)

Specifications are subject to change without prior notice.

In case of any questions, please contact our sales representative via phone, e-mail or our website.

<https://www.ritter-medical.de/en/contact>