


True to CAPP's commitment to quality, our new petri dishes are manufactured to meet industry-leading quality measures. Manufactured in a contaminant-free cleanroom and from premium-quality polystyrene, CAPP Petri Dishes offer truly excellent performance.

The industrial design of each petri dish allows for maximum heat transfer between the ambient air in the incubator and the bacteria on the culture. This creates an optimal thermal zone for growth at the set temperature. With superior optical properties, tasks like observing growth patterns, hemolysis and antibiotic sensitivity, among others, are quick and accurate.

Each batch is supplied in packs of 10 or 20 petri dishes. Users can open the right number of plates for their culture media preparation without wasting media or contaminating unused plates. These petri dishes can be ordered in vented and non-vented formats, both of which optimise airflow for different culture conditions. All CAPP Petri Dishes have a standard 90mm diameter that provides ample room for streak-plates or spread-plates. For space optimisation, users can stack CAPP Petri Dishes in an incubator or refrigerator for long-term storage.

### Features:

- Excellent optical properties
- Space-saving stackable design
- Premium-quality polystyrene construction
- Contaminant-free petri dishes
- 90mm wide surface area for plating
- Simple packaging with traceable lot numbers



## Ordering Information

Cat. No.	Description
5061000C	CAPP Petri Dishes, Ø91.5mm, H:14.35mm, w/o vents, non-sterile, 48x10 pcs.
5061001C	CAPP Petri Dishes, Ø91.5mm, H:14.35mm, w/o vents, sterile, 48x10 pcs.
5061010C	CAPP Petri Dishes, Ø91.5mm, H:14.55mm, with 3 vents, non-sterile, 48x10 pcs.
5061011C	CAPP Petri Dishes, Ø91.5mm, H:14.55mm, with 3 vents, sterile, 48x10 pcs.
5062000C	CAPP Petri Dishes, Ø91.5mm, H:14.35mm, w/o vents, non-sterile, 24x20 pcs.
5062001C	CAPP Petri Dishes, Ø91.5mm, H:14.35mm, w/o vents, sterile, 24x20 pcs.
5062010C	CAPP Petri Dishes, Ø91.5mm, H:14.55mm, with 3 vents, non-sterile, 24x20 pcs.
5062011C	CAPP Petri Dishes, Ø91.5mm, H:14.55mm, with 3 vents, sterile, 24x20 pcs.

